



ProLyte®

The World's Most Popular Electrolyte Analyzer

 [diamonddiagnostics.com/prolyte](https://www.diamonddiagnostics.com/prolyte)

ProLyte®

The World's Most Popular Electrolyte Analyzer



Sample Size

Whole Blood, Serum (95 µl)
Plasma, QCs (95 µl)
Urine (180 µl)

Languages

English, 中文, Español, Français,
Deutsch, Bahasa Indonesia, Italiano,
日本語, 한국어, Português, Polskie,
Русский, Türkçe

Analysis Time

45 sec (240/hour)

User Input

4" Touch Screen
External Keyboard
Remote Access

Number of Parameters†

Na⁺ K⁺ Cl⁻ Li⁺
4 Total
3 Simultaneously

Connectivity

LAN, Wifi,
4 USB ports
RS232, LIS

Data Storage

20,000 Patient Results
10,000+ QC Results per Level

Calibration

2 Point every 4 hours
1 Point every sample
On Demand

Support

Online Service Interface

Parameter	Range	Reproducibility ¹	Resolution
Na ⁺	45-205 mmol/L	CV ≤ 1% (120-160mmol/L)	0.01 mmol/L
K ⁺	1.5-15 mmol/L	CV ≤ 1.5% (2.5-6 mmol/L)	0.01 mmol/L
Cl ⁻	45-205 mmol/L	CV ≤ 1% (85-130 mmol/L)	0.01 mmol/L
Li ⁺	0.15-5.00 mmol/L	SD ≤ 0.03 mmol/L (0.4-1.3 mmol/L)	0.001 mmol/L
Urine ² Na ⁺	30-1020 mmol/L	CV ≤ 5% (100-250 mmol/L)	0.01 mmol/L
K ⁺	20-505 mmol/L ³	CV ≤ 5% (10-60 mmol/L)	0.01 mmol/L
Cl ⁻	25-505 mmol/L	CV ≤ 5% (100-250 mmol/L)	0.01 mmol/L

1- Typical Within Run (n=30) Blood, Serum, Plasma
2- Lithium is not typically measured in urine samples
3- (60-120) requires additional dilution

†Calcium available for certain markets. Not available in USA.

Power

100-240V, 50W max
50-60 Hz, 1.6A max

Size & Weight

13.2 x 12.4 x 12 in, 14 lbs
335 x 315 x 295 mm, 6 kgs

Ambient Conditions

Humidity: < 85%
Room Temperature: 60-90°F (15-32°C)



DIAMOND
Smart Lab Solutions

333 Fiske Street, Holliston, MA 01746 U.S.A.
www.diamonddiagnostics.com
main@diamonddiagnostics.com
+1 508-429-0450

SP023-0005F REV29 01/22/25